



## BS-Seq sample preparation guidelines

### Genomic DNA for WGBS-Seq and RRBS-Seq. Input Requirements.

IGATech offers nucleic acids extraction service and can set up a dedicated extraction workflow for a specific substrate. Please enquire.

#### DNA Quantity

Submit 1 µg of DNA per sample at a minimum concentration of 20 ng/µL.

Please note that fluorimetry-based quantification (*e.g.* Qubit, plate-reader) assays are more accurate than absorbance-based methods (*e.g.* Nanodrop), which might overestimate the quantity.

#### DNA Purity

For WGBS-seq experiments, resuspend DNA in low-EDTA TE buffer or Nuclease-free Water.

In case of RRBS-seq application, DNA must be submitted in Nuclease-free Water.

DNA samples must be free of contaminating proteins, RNA, organic solvents (including phenol and ethanol) and salts.

The A260:A280 and A260:A230 ratios for DNA samples should be > 1.8. Use of DNA with lower ratios may result in low amplification yield.



## DNA Integrity

Prepare the DNA following your favorite extraction method, even if it is strongly recommended to use commercial column-based protocols.

Control DNA on 0.8% agarose gel to check for integrity. Degraded gDNA may affect the quality of the final libraries, leading to over-fragmentation of DNA and insufficient libraries complexity.

**For challenging samples (low inputs or FFPE specimens) please contact us to determine the feasibility of the processing.**

## Shipping

**For batches of <24 samples**, send samples in 1.5 mL or 2 mL Eppendorf tubes sealed with parafilm (***0,5 mL and 0.2 mL tubes as well as strips will be not accepted***). The tubes must have, on the vial top, a clear and permanent sign (or a thin label) with a **progressive number** corresponding to information specified in the **Sample Spreadsheet**.

**For batches of 24 or more samples**, send samples in a skirted 96-wells plate, sealed with adhesive/heat-sealed aluminum foil. Each plate must be labeled with a plate identifier indicated in the **Sample Spreadsheet**.

Ship samples in a cold pack (*e.g.* Blue ice). Organize the packaging in a way to avoid damaging and dispersion during the shipment: place tubes and plates in smaller plastic bags or boxes.

It is **MANDATORY** to send us the compiled **Sample Spreadsheet**, both with the shipped parcel and via e-mail. To be able to properly track and safeguard your samples send us the **Tracking Number** via e-mail.



### Summary of best practice and quality requirements for NGS

260/280 > 1.8 and 260/230 > 1.8

Quantify by dsDNA-specific fluorimetry (Qubit/fluorimeter)

Avoid repeated freeze-thaw cycles (arrange aliquots)

Do not expose to high temperatures (>65 C°) especially for HMW DNA

Store the DNA in stabilizing buffer (Tris-HCl pH 8.0-8.5)

Avoid EDTA (conc. > 0.1 mM)

DNA must be RNA-free (RNase treatment is strongly recommended)

DNA must be free of phenol and polyphenols