






GUIDANCE ON SAMPLE PLATES DELIVERY AND SHIPMENT



Projects with fewer than 24 samples




Labelled **1.5/2.0 ml tubes** with safe-lock lids
or parafilm



0.2- & 0.5-mL PCR strip or tubes

Projects with fewer than 24 samples



- Label tubes simply with progressive numbers using a waterproof marker pen
- Pay attention that the tube is perfectly closed (i.e., safe-lock lids, parafilm)



- Avoid hardly visible marker pen (better black and blue, rather than orange or pink)
- Avoid extra information in the label (e.g., sample name, date, extraction kit)

Projects with fewer than 24 samples

Protection from damage by use of right size boxes, bubble wrap or similar.







Projects with more than 24 samples

Please submit your samples in a **well-sealed and labelled 96-well plate**

Plates must be one of the following or similar in their specification and sizes:


- *MicroAmp Optical 96-Well Reaction Plate' (Thermo Fisher™) REF N8010560*
- *Eppendorf 96-Well twin.tec™ PCR Plates (Eppendorf™) Order no: 0030 128.672, Cat. No: 951020460*
- *Sarstedt PCR plate half skirt, 96 well, transparent, High-Profile, 200 µl, PCR Performance Tested, PP; Order no: 72.1979*




Projects with more than 24 samples

Please submit your samples in a **well-sealed and labelled 96-well plate**

PER-COLUMN
ORDER





PER-ROW
ORDER



Projects with more than 24 samples

Please place the plate in the insert from a tip box to prevent disruption and/or distortion



Projects with more than 24 samples

Please submit your samples in a **well-sealed and labelled 96-well plate**



PCR Plate Sealing Mat



PCR strip caps, closing perfectly each well.
Protect your plate wrapping it in **parafilm**

WITHOUT HEATSEALER





PCR Sealing Film



Aluminum Foil

Projects with more than 24 samples

In absence of a heatsealer, ensure the seal is completely attached along the four sides and around each well to prevent contamination of and between samples as well as evaporation.



Sample Spreadsheet preparation

Fill-in all the contact information
In particular the quotation number

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
1	Quote*				QUO#####											
2	Reference Person*				"the PI name"											
3	Contact Person*				"the sender name"											
4	Contact e-mail* (recipient of work in progress reports)				"the mail of the sender"											
5	University/Department				please also put your organization here to help sample identification											
6	Company/Institute				ase also put your organization here to help sample identification											
7	Telephone				not strictly necessary											
9																
10	Samples definition														Bioinformatics & Computational	
11	Human-derived samples must be anonymized. Do not provide sensitive data.														Meta data fields for comparison/statistical experiments -- categorical data	
12	For samples storage and return policies, please refer to our Terms and Conditions. https://igatechtechnology.com/igatech/documents/														Mostly required for Metagenomics, RNA-Seq, ChIP-Seq STANDARD A	
13	Please do not paste sample names with space or special characters, use max 12 characters														Please inquire in case of other applications.	
14	* mandatory field - DO NOT USE SPECIAL CHARACTERS ** mandatory for multi-plates ONLY															
15	Do not use this column														Do not use this column	
16	do not compile	Customer label* (max 12 characters)	Plate identifier**	Well position** in per-column order	Sample type*	SPECIES (latin binomial)*	VOLUME (ul) *	CONCENTRATION (ng/ul) *	AMOUNT (ug)	A260 A280	A260 A230	R.I.N. D.I.N.	Storage Buffer or Extraction kit*	IDXXXX	Condition	Experimental plan
17	1	CS1	P1	A1	DNA	Microbial metasamp	20	1.3	0.03	1.56	0.45		TE buffer	IDXXXX_1-CS1-P1-A1	Group1	Group1_vs_Group2
18	2	CS2	P1	A2	DNA	Microbial metasamp	20	3.3	0.07	1.34	0.19		TE buffer	IDXXXX_2-CS2-P1-A2		Group1_vs_Group3
19	3	CS3	P1	A3	DNA	Microbial metasamp	20	1.3	0.03	1.97	0.23		TE buffer	IDXXXX_3-CS3-P1-A3		Group2_vs_Group3
20	4	CT1	P1	A4	DNA	Microbial metasamp	20	14.9	0.30	1.38	0.08		TE buffer	IDXXXX_4-CT1-P1-A4	Group2	
21	5	CT2	P1	A5	DNA	Microbial metasamp	20	7.3	0.15	1.58	0.05		TE buffer	IDXXXX_5-CT2-P1-A5		
22	6	CT3	P1	A6	DNA	Microbial metasamp	20	4.8	0.10	1.88	0.05		TE buffer	IDXXXX_6-CT3-P1-A6		
23	7	WS1	P1	A7	DNA	Microbial metasamp	20	1.3	0.03	0.87	0.07		TE buffer	IDXXXX_7-WS1-P1-A7	Group3	
24	8	WS2	P1	A8	DNA	Microbial metasamp	20	3	0.06	2.22	0.12		TE buffer	IDXXXX_8-WS2-P1-A8		
25	9	WS3	P1	A9	DNA	Microbial metasamp	20	4.3	0.09	1.34	0.39		TE buffer	IDXXXX_9-WS3-P1-A9		

Avoid spaces and special characters in the
sample name: only use numbers and
letter separated by dash "-"

Box preparation guidelines


Prevent damage of individual wells by use of protective wrapping and place the plate in a bag within box.
Ensure that sample tubes and plates are protected from damage by use of bubble wrap or similar.



Box preparation guidelines


Specific precautions for **RNA samples**

Send RNA samples in **dry-ice** or **lyophilized** with GenTegraRNA or GENEWIZ RNA Stabilization Tubes.



Box preparation guidelines

Ensure that sample tubes and plates are well fixed in shipping box, enclose the **Sample Spreadsheet** and close your box. **You are now ready for delivery**



Send the **Sample Spreadsheet** file via email to either:

- your IGATech contact
- info@igatechnology.com

Thank You

