



## Metabarcoding sample preparation guidelines

IGATech offers nucleic acids extraction service, and we can set up a dedicated extraction workflow for your specific substrate. Please enquire.

### Amplicon specifications

In case you are sending amplicons, **specify in the Sample Spreadsheet whether Illumina overhang adapter sequences have been appended to locus-specific sequences** (primers):

Forward overhang: 5' TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG-[*locus-specific sequence*]

Reverse overhang: 5' GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG-[*locus-specific sequence*]

More details can be found on page 3 of

[https://support.illumina.com/documents/documentation/chemistry\\_documentation/16s/16s-metagenomic-library-prep-guide-15044223-b.pdf](https://support.illumina.com/documents/documentation/chemistry_documentation/16s/16s-metagenomic-library-prep-guide-15044223-b.pdf)

**In case Illumina adapters haven't been used please make it clear in the quotation request since the per sample price is different.**

Make sure to **specify amplicon length** and **protocol used for purification**. In case you're sending a mix of amplicons specify also if their sequences are similar or highly variable.

We suggest sending purified amplicon in water or 10mM Tris-HCl pH 8.5 at a minimum concentration of 10 ng/ $\mu$ L, min volume 25 $\mu$ L. For amplicons longer than 500 bp send at least 500 ng in min volume 25 $\mu$ L.

See the **shipping information** below.

### DNA specifications for 16S - 18S – ITS

We suggest sending 250 ng of DNA at a minimum concentration of 10 ng/ $\mu$ L, min volume 25 $\mu$ L. Send DNA samples re-suspended in water or 10mM Tris-HCl pH 8.5.

Sample acceptance is based on the amount and concentration measured by fluorimetry-based quantification (*e.g.*, Qubit, plate-reader). Keep in mind that absorbance-based methods (*e.g.*, Nanodrop) might largely **overestimate** the DNA quantity.



The A260:A280 and A260:A230 ratios for DNA samples should be  $> 1.8$ . Use of DNA with lower ratios may result in low amplification yield.

We understand that from some substrates it can be hard to obtain high yields of DNA. We also accept samples below the recommended quantity/purity. It is recommended, when possible, that customers run a test PCR to ensure that the samples are reaction-permissive, as the DNA of “dirty” substrates (soil, sludge, fecal, etc.) may contain PCR-inhibitor leftovers which can hamper amplification reaction.

## Shipping

**For batches of <24 samples**, send samples in 1.5 mL or 2 mL Eppendorf tubes sealed with parafilm (**0,5 mL and 0.2 mL tubes as well as strips will be not accepted**). The tubes must have on the vial top, a clear and permanent sign (or a thin label) with a **progressive number** corresponding to information specified in the **Sample Spreadsheet**.

**For batches of >24 samples**, send samples in a skirted 96-wells plate, sealed with adhesive/heat-sealed aluminum foil. Each plate must be labeled with a plate identifier indicated on the **Sample Spreadsheet**.

Ship samples in a cold pack (e.g., **Blue ice**). Organize the packaging in a way to avoid damaging and dispersion during the shipment: place tubes and plates in smaller plastic bags or boxes. Put a separator between stacked plates to avoid perforations of adhesive foils and leakage.

Significant differences in nucleic acid concentration of samples can occur if samples are of different origin or extracted with different methods. This may cause different efficiencies during library construction. Please, organize plates by creating homogeneous groups and report all information in the **Sample Spreadsheet**.

It is **MANDATORY to send us the compiled Sample Spreadsheet**, both with the shipped parcel and via e-mail. To be able to properly track and safeguard your samples send us the Tracking Number via e-mail.

**Customer is free to add extra columns to Sample Spread Sheet to include all available metadata of samples, if the analysis service has been requested.** This information will be used to run several statistical tests as well as labeling/coloring of outputs in our standard analysis.



***IMPORTANT: We remind you that human-derived samples must be anonymized. Therefore, we cannot accept samples that come along with personal identification data (name and surname, fiscal code, etc.). Supplementary data related to the study such as prognosis, biometrics values, age, sex, and other information not directly associated to an individual can be provided with no limitation (please use the Sample Spreadsheet available in the Documents section).***

***If human samples are provided in the form of tissue or body fluid, please fill out the [Human samples clearance form](#) and return a signed copy.***